

## PROCEDURE FOR SOLVENT CLEANING OF SILICON WAFERS

1. Put the wafer on top of the spin coater chuck located in the fume hood-wet bench.



2. Check recipe conditions for:
  - a. Spin speed range between 2000 and 3000rpm
  - b. Total time of 8 mins
3. For 1 min, pour Acetone from the squeezable bottle on the rotating wafer and use foam-cotton swab to remove gently any residues. Do not stop pouring acetone during this lapse of time.
4. For 1 min, pour Methanol from the squeezable bottle on the rotating wafer and use foam-cotton swab to remove gently any residues. Do not stop pouring acetone during this lapse of time.
5. For 1 min, pour Isopropyl Alcohol from the squeezable bottle on the rotating wafer and use foam-cotton swab to remove gently any residues. Do not stop pouring acetone during this lapse of time.
6. Rinse wafer with DI water during 1 min.
7. Use pressurized nitrogen to blow dry for 1 min.
8. Turn off the spinner vacuum and proceed to remove the wafer from the chuck using tweezers to avoid touching the clean surface.